


SQUARE ABUTMENTS



bti.[®]

Human
Technology



KNOW THE MANY ADVANTAGES OF THE NEW

SQUARE ABUTMENTS

- Enables you to create a **BIO emergence profile** suitable for the needs of each patient.
- **Simplifies the manufacture of abutments** in terms of time and experience.
- **Hermetic implant-prosthesis seal.**

CONVENTIONAL AND DIGITAL TECHNIQUES

Thanks to their versatility, they can be used with or without a castable sleeve, both for the production of conventional prostheses and for CAD/CAM systems (they are used as machined bases, on which Zr, e-max, Ti, Cr-Co and resin copings can be modelled).

CLINICAL INDICATIONS

1. Single and multiple direct to implant prosthesis.
2. Immediate loading of single implants in aesthetic area.



SQUARE ABUTMENTS: EXCELLENT AESTHETICS AND GREAT VERSATILITY TO OPTIMISE EVERY CASE

Innovative range of square abutments for direct-to-implant prostheses, featuring a BIO profile for long-term maintenance of the gingival mucosa, specific surface treatment that increases the biological seal and the possibility of changing the angulation at screw emergence.

► MATERIAL

Titanium (Ti) with Ti-Golden® surface treatment.

► BIO PROFILE

BIO prosthetic profile for maintaining the gingival vascularisation and minimising the risk of recession in cases of fine biotypes.

► CASTING SLEEVES

- For conventional lab processes: ceramic injection and casting.
- For narrow and universal platforms.



► AVAILABLE PLATFORMS: 3.0, narrow and universal.

► VARIOUS HEIGHTS

The square abutments are available in various gingival heights: 0,5 - 1 - 2 and 3 mm.



► SPECIFIC SCREW (INCLUDED)

Ensure preloading.

► IMPROVEMENTS AND INNOVATION

- Slotted chimney for improved retention.
- Two customisation heights: 4 and 6,5 mm for better adaptation to the prosthetic space.
- Side window to allow angulations up to 15°.

► ANTI-ROTATIONAL DESIGN

Facilitates the prosthetic work in CAD/ CAM.

► CONNECTION TO THE IMPLANT


Non rotational

Rotational



Single and multiple cement-retained prostheses

Multiple screw-retained prosthesis.



Products manufactured by BTI comply with the current Medical Device Directive 93/42/EEC, the Medical Device Regulation 2017/745, the FDA Quality System Regulation 21 CFR Part 820, the Canadian Medical Device Regulation SOR/98-282 and the Medical Device Quality Management System ISO 13485. The Notified Body TUV SUD, No. 0123, has assessed the conformity of these products.



BTI Comercial

San Antonio, 15 - 5º
01005 Vitoria-Gasteiz
(Álava) · SPAIN
Tel: +34 945 140 024
pedidos@bticomercial.com

B.T.I.

Biotechnology Institute S.L.
Parque Tecnológico de Álava
Leonardo da Vinci ,14
01510 Miñano (Álava) SPAIN
bti.implantes@bti-implant.es

Subsidiaries

GERMANY

Tel: +49 7231 428060
info@bti-implant.de

PORTUGAL

Tel: (351) 22 120 1373
bti.portugal@bticomercial.com

FRANCE

Tel: +33 (5) 54530209
serviceclient@bti-implant.fr

UK

Tel: +44 02039 661873
customerservice@bti-implant.co.uk

ITALY

Tel: +39 0270605067
info@bti-implant.it

USA

Office: 866 646 4067
bti@bti-implant.us

MEXICO

Tel: (52) 55 52502964
bti.mexico@bti-implant.com

www.bti-biotechnologyinstitute.com

NOTE: consult your distributor for the availability of the product in the different markets.